





Birmingham Quality


Serum creatinine Median B is +0.7 and C is 2.9 for Compensated Kinetic Jaffe Abbott reagents [10AB]




Serum creatinine Median B is +1.8 and C is 4.9 for Compensated Kinetic Jaffe Roche Cobas/Modular [10BO]




Serum creatinine Median B is -3.9 and C is 3.7 for Compensated Kinetic Jaffe Olympus reagents [10OL]




Serum creatinine Median B is -5.2 and C is 4.3 for Compensated Kinetic Jaffe Siemens ADVIA [10TE]




Serum creatinine Median B is +0.6 and C is 3.8 for Dry slide OCD (J&J) Gen 80+ [1JJ2]



Serum creatinine Median B is +0.8 and C is 1.6 for Enzymatic Abbott reagents [9AB]



Serum creatinine Median B is +3.1 and C is 2.5 for Enzymatic Roche Cobas/Modular [9BO]



Serum creatinine Median B is -8.5 and C is 3.9 for Enzymatic Siemens ADVIA [9TE]

